

INTENDED USE


The **Carica papaya Papain ELISA Kit** is a sandwich ELISA for the detection and/or quantification of papain in solution. The kit is designed to measure trace residue of papain in vaccine, therapeutic antibodies, drugs or other biological samples, regardless of its enzyme activity; for research use only (RUO), not for diagnosis, cure or prevention of the disease.

GENERAL INFORMATION



Papain, also known as papaya proteinase I, is a cysteine protease (EC 3.4.22.2) enzyme present in papaya (*Carica papaya*) and mountain papaya (*Vasconcellea cundinara maricensis*). The papain

family of enzymes shows a wide variety of activities, including endopeptidases, aminopeptidases, dipeptidyl peptidases and enzymes with both exo- and endo-peptidase activity. Papain related enzymes are widespread, found in baculovirus (insect), eubacteria, yeast, and practically all protozoa, plants and mammals. The proteins are typically lysosomal or secreted, and proteolytic cleavage of the propeptide is required for enzyme activation. Papain-like cysteine proteinases are essentially synthesised as inactive proenzymes (zymogens) with N-terminal propeptide regions. The activation process of these enzymes includes the removal of propeptide regions, which serve a variety of functions in vivo and in vitro. Amino acid residues within the pro-region mediate their membrane association, and play a role in the transport of the proenzyme to lysosomes.




The papain precursor protein contains 345 amino acid residues, and consists of a signal sequence (1-18), a propeptide (19-133) and the mature peptide (134-345). The amino acid numbers are based on the mature

peptide. The protein is stabilised by three disulfide bridges. Papain is a relatively heat-resistant enzyme, with an optimal temperature range of 60 and 70 °C. Papain breaks down tough meat fibres, and has been used for thousands of years to tenderise meat. Meat tenderisers in powder form with papain as an active component are widely sold. Papain can be used to dissociate cells in the first step of cell culture preparations. Papain is added to some toothpastes and mint sweets as a tooth whitener. Topical drug ointments (Accuzyme, Allafil, Allanzyme, Ethezyme, Gladase, Kovia, Panafil, Pap Urea, and Ziox) containing papain are used to remove dead or contaminated tissue in acute and chronic lesions, such as diabetic ulcers, pressure ulcers, varicose ulcers, and traumatic infected wounds. Other products are marketed under the names of the active ingredients, for instance, papain-urea ointment. In 2008 the FDA announced its intention to take action against these products because it had received reports of serious adverse events in patients using products containing papain.

People allergic to latex can also be allergic to papaya, the source of papain, implying that people with latex sensitivity may be at increased risk of suffering an adverse reaction to a topical papain drug product. Therefore, sensitive immunoassay for papain are required to detect the presence of the enzyme even after inactivation.

An antibody digested by papain yields three fragments: two 50 kDa Fab fragments and one 50kDa Fc fragment. The papain-digested antibody is unable to promote agglutination, precipitation, opsonization, and lysis. Papain is used for fragmentation of many therapeutic antibodies (CroFab, DigiFab, DigiBind, ReoPro, Lucentis and Cimzia). Papain is routinely used for despeciation of anti-venom produced in horse, goat, and sheep.

PRINCIPLE OF THE TEST



The Human Beta Defensin 2 ELISA kit is based on the binding of Human Beta Defensin 2 in samples to two antibodies, one immobilized on the microtiter wells, and the other conjugated to biotin, which then binds to a streptavidin horseradish peroxidase (HRP) conjugate. After a washing step, chromogenic substrate is added and color is developed by the enzymatic reaction of HRP on the TMB substrate, which is directly proportional

to the amount of Papain present in the sample. Stopping Solution is added to terminate the reaction, and absorbance at 450nm is then measured using an ELISA microtiter well reader. The concentration of Papain in samples is calculated from a standard curve of purified recombinant human BD-2 of designated concentration.

STORAGE AND STABILITY

The microtiter well plate and all other reagents, if unopened, are stable at 2-8°C until the expiration date printed on the label. Stabilities of the working solutions are indicated under Reagent Preparation.

KIT CONTENTS

Ready For Use: Store as indicated on labels.

Component	Part #	Amt	Contents
Anti-Papain Coated Strip Plate	800-161	8-well strips (12)	Coated with purified anti-papain antibodies. Return unused strips to the pouch with desiccant; re-seal and store refrigerated.
TMB Substrate	80091	12 ml	Chromogenic substrate for HRP containing TMB and peroxide.
Stop Solution	80101	12 ml	1% sulfuric acid.

Materials Required But Not Provided:

- Pipettors and pipettes that deliver 100ul and 1-10ml. A multi-channel pipettor is recommended.
- Disposable glass or plastic 5-15ml tubes for diluting samples and Antibody HRP Concentrate.
- Graduated cylinder to dilute Wash Concentrate; 0.2 to 1L.
- Stock bottle to store diluted Wash Solution; 200ml to 1L.
- Distilled or deionized water to dilute reagent concentrates.
- Microwell plate reader at 450 nm wavelength.

To Be Reconstituted or Diluted: Store as indicated.

Component	Instructions for Use																												
Papain Standard Part No. 800-162	Three (3) vials, each containing Papain lyophilized in buffer with proteins and antimicrobial. Refrigerate lyophilized vials until used or kit lot expires.																												
Reconstitute 1 vial with the volume of Working Sample Diluent indicated on the Standard label to provide a 100 ng/ml Top Standard. Prepare 2-fold dilutions, as follows; sufficient for one entire curve:																													
	<table border="1"> <thead> <tr> <th>Standard</th> <th>+ Diluent</th> <th>=</th> <th>Final Conc</th> </tr> </thead> <tbody> <tr> <td>Reconstituted Standard</td> <td>None</td> <td></td> <td>100 ng/ml</td> </tr> <tr> <td>225 ul of</td> <td>100 ng/ml</td> <td>225ul</td> <td>50 ng/ml</td> </tr> <tr> <td>225 ul of</td> <td>50 ng/ml</td> <td>225ul</td> <td>25 ng/ml</td> </tr> <tr> <td>225 ul of</td> <td>25 ng/ml</td> <td>225ul</td> <td>12.5 ng/ml</td> </tr> <tr> <td>225 ul of</td> <td>12.5 ng/ml</td> <td>225ul</td> <td>6.25 ng/ml</td> </tr> <tr> <td>225 ul of</td> <td>6.25 ng/ml</td> <td>225ul</td> <td>3.125 ng/ml</td> </tr> </tbody> </table>	Standard	+ Diluent	=	Final Conc	Reconstituted Standard	None		100 ng/ml	225 ul of	100 ng/ml	225ul	50 ng/ml	225 ul of	50 ng/ml	225ul	25 ng/ml	225 ul of	25 ng/ml	225ul	12.5 ng/ml	225 ul of	12.5 ng/ml	225ul	6.25 ng/ml	225 ul of	6.25 ng/ml	225ul	3.125 ng/ml
Standard	+ Diluent	=	Final Conc																										
Reconstituted Standard	None		100 ng/ml																										
225 ul of	100 ng/ml	225ul	50 ng/ml																										
225 ul of	50 ng/ml	225ul	25 ng/ml																										
225 ul of	25 ng/ml	225ul	12.5 ng/ml																										
225 ul of	12.5 ng/ml	225ul	6.25 ng/ml																										
225 ul of	6.25 ng/ml	225ul	3.125 ng/ml																										
Stable for up to 4 weeks at 2-8°C.																													
Sample Diluent Concentrate (20x) Cat. No. SD-20T, 10ml	Dilute the entire volume, 10ml + 190ml with distilled or deionized water into a clean stock bottle. Label as Working Sample Diluent and store at 2-8°C until the kit lot expires or is used up.																												
Wash Solution Concentrate (100x) Cat. No. WB-100, 10ml	Dilute the entire volume, 10ml, to 1L with distilled or deionized water into a clean stock bottle. Label as Working WashSolution and store at ambient temperature until kit is used entirely.																												
Anti-Papain-HRP Conjugate Concentrate (100x) Part No. 800-164, 0.15ml	Peroxidase conjugated anti-papain in buffer with protein, detergents and BND as stabilizers. Dilute fresh as needed; 10ul of concentrate to 1ml of Working Sample Diluent is sufficient for 1 8-well strip. Use within the working day and discard. Return concentrate to 2-8°C storage.																												

ASSAY DESIGN AND SET-UP

Sample Collection and Handling

Culture medium, bioprocessing preparations, serum and other biological fluids may be used as samples with proper dilution to avoid solution matrix interference (See Limits of the Assay, page 6). For all samples, clarify by centrifugation and/or filtration prior to dilution in Sample Diluent. If samples will not be assayed immediately, store refrigerated for up to a few weeks, or frozen for long-term storage.

Assay Validation

Validate the performance of the sample antigen and matrix in the assay system for recovery and parallelism (see Limits of the Assay, page 6), as follows:

Recovery – a measure of the interference of the sample matrix (diluent effect) in providing accurate quantitation of the sample Humira relative to the Standard curve.

Prepare and run a series of dilutions of the sample antigen (concentrations that will fall within the Standard range) in Working Sample Diluent to determine the dilutions that give consistent and accurate quantitation. For most buffer solutions a minimum 5-fold sample dilution is usually sufficient. Serum and plasma require at least a 10-fold dilution to obtain consistent quantitation or complete antigen recovery.

Parallelism – dilutions of the sample should read equivalent values from the top and bottom of the Standard curve to provide good assay precision.

Prepare a dilution series of the sample antigen that gives complete recovery and falls within the full range of the Standard curve. Sample readings from the upper and lower regions of the curve should differ by less than 25%.

Plate Set-up

Bring all reagents to room temperature (18-30° C) equilibration (at least 30 minutes).

- Determine the number of wells for the assay run. Duplicates are recommended, including 10 Standard wells and 2 wells for each sample and control to be assayed.
- Remove the appropriate number of microwell strips from the pouch and return unused strips to the pouch. Reseal the pouch and store refrigerated.
- Prior to sample addition, add 200-300ul Working Wash Solution to each well and let stand for about 5 minutes.
- Aspirate or dump the liquid and pat the plate dry on a paper towel.

Assay Procedure

ALL STEPS ARE PERFORMED AT ROOM TEMPERATURE. After each reagent addition, gently tap the plate to mix the well contents prior to beginning incubation.

- 1st Incubation [100ul – 60 min; 4 washes]**
 - Add 100ul of calibrators, samples and controls each to pre-determined wells.
 - Tap the plate gently to mix reagents and incubate for 60 minutes.
 - Wash wells 4 times and pat dry on fresh paper towels. As an alternative, an automatic plate washer may be used. Improper washes may lead to falsely elevated signals and poor reproducibility.
- 2nd Incubation [100ul – 30 min; 5 washes]**
 - Add 100ul of diluted Anti-papain HRP to each well.
 - Incubate for 30 minutes.
 - Wash wells 5 times as in step 2.
- Substrate Incubation [100ul – 15 min]**
 - Add 100ul TMB Substrate to each well. The liquid in the wells will begin to turn blue.
 - Incubate for 15 minutes in the dark, e.g., place in a drawer or closet.

Note: If your microplate reader does not register optical density (OD) above 2.0, incubate for less time, or read OD at 405-410 nm (results are valid).
- Stop Step [Stop: 100ul]**
 - Add 100ul of Stop Solution to each well.
 - Tap gently to mix. The enzyme reaction will stop; liquid in the wells will turn yellow.
- Absorbance Reading**
 - Use any commercially available microplate reader capable of reading at 450nm wavelength. Use a program suitable for obtaining OD readings, and data calculations if available. Read absorbance of the entire plate at 450nm within 30 minutes after Stop Solution.

CALCULATION OF RESULTS


The results may be calculated using any immunoassay software package. The four-parameter curve-fit is recommended. If software is not available, Papain concentrations may be determined as follows:

1. Calculate the mean OD of duplicate samples.
2. On graph paper plot the mean OD of the standards (y-axis) against the concentration (ng/ml) of Papain (x-axis). Draw the best fit curve through these points to construct the standard curve. A point-to-point construction is most common and reliable.
3. The Papain concentrations in unknown samples and controls can be determined by interpolation from the standard curve.
4. Multiply the values obtained for the samples by the dilution factor.
5. Samples producing signals higher than the 200 ng/ml standard should be further diluted and re-assayed.

Typical Results:

Wells	Standards, Control & Samples	A450	ng/ml
A1, A2	Negative Diluent Control	0.04	0
B1, B2	3.125 ng/ml Standard	0.30	3.125
C1, C2	6.25 ng/ml Standard	0.55	6.25
D1, D2	12.5 ng/ml Standard	0.89	12.5
E1, E2	25 ng/ml Standard	1.36	25
F1, F2	50 ng/ml Standard	1.98	50
G1, G2	100 ng/ml Standard	2.97	100
H1, H2	Sample [Diluted 1:10]	1.55	31

Calculated: 10-fold dilution x 31 ng/ml = **310** ng/ml in sample



PERFORMANCE CHARACTERISTICS

Specificity

The antibodies used in this kit have been affinity-purified using a papain immunosorbent and have been shown by immunoelectrophoresis and ELISA to react specifically with papain.

Precision

Samples containing low, medium and high concentrations of Papain as duplicates in multiple assays (n=5) to obtain between-assay reproducibility. Coefficients of variation were calculated for the concentrations using a point-to-point curve-fitting program.

Papain concentrations were measured with good between-assay (4.3 to 9.2 %CV) reproducibility.

Sample	Papain ng/ml	Inter-assay % CV
High Concentration	34.1	9.2
Medium Concentration	18.4	4.3
Low Concentration	8.6	6.8

References: Cohen LW (1986) Gene 48, 219-227; Dretnh J (1968) Nature 218, 929-932; Rawligns ND (1994) Methd. Enzymol. 244, 461-486; Sebt SM (1987) Biochem. J. 290, 205-218; Yamamoto Y (2002) *Curr Protein Pept Sci* 3 (2): 231-238; Shuren J (2008) "Topical Drug Products Containing Papain; Enforcement Action Dates" (PDF). United States Food and Drug Administration, Department of Health and Human Services.

LIMITS OF THE ASSAY

1. Papain that is incomplete in sequence (truncated) or is aggregated and/or associated with other biomolecules may not produce dilution curves **parallel** with the Standard curve. For cases of non-parallelism, it may be useful to establish an alternative Standard curve using the altered papain preparation.
2. Papain in an activated buffer may require de-activation prior to its use as a sample in the assay, or proteolysis of the coating antibodies could produce falsely low signals. Generally, inactivation by removing sulfhydryl activators (e.g., by dialysis) is sufficient, but must be verified by the operator.

QUALITY CONTROL

Reagents Accurate and reproducible assay results rely on proper storage, handling and control of reagent and sample temperature. Store all reagents as indicated, and warm to room temperature only those to be used in the assay. Shelf-life of the critical reagents and samples will diminish with extended exposure to non-refrigeration, resulting in inaccurate assay results. All solutions should be clear. Cloudiness or particulates are indications of reagent contamination or instability and may interfere with proper performance of the assay. Do not use.

Standard Curve The signal generated by the standards should be continuously increasing in OD from the lowest Standard to the highest Standard, with a difference greater than 1.2 OD. Non-uniform or low signals may indicate problems with technique, protocol directions and/or reagent preparation, use or stability. Do not rely on results generated from an assay with these issues.

Technique Accurate and reproducible assay results rely on good lab technique regarding pipetting, plate washing and handling of samples and reagents.

Equipment Precision of results relies on uniform and effective washing techniques; an automatic washer may be used. ELISA reader and pipettes should be properly calibrated.

PRECAUTIONS AND SAFETY INSTRUCTIONS

Standards, Sample Diluent, and Antibody HRP contain bromonitrodioxane (BND: 0.05%, w/v). Stop Solution contains 1% sulfuric acid. Follow good laboratory practices, and avoid ingestion or contact of any reagent with skin, eyes or mucous membranes. All reagents may be disposed of down a drain with copious amounts of water. MSDS for TMB, sulfuric acid and BND can be requested or obtained from the ADI website: http://4adi.com/commerce/info/showpage.jsp?page_id=1060&category_id=2430&visit=10

Papain ELISA Kit

Cat. 800-160-CPP, 96 Tests

For Quantitation of Papain in Solution

For research use only (RUO), not for diagnostic or therapeutic use.



**ALPHA DIAGNOSTIC
INTERNATIONAL**

6203 Woodlake Center Drive • San Antonio • Texas 78244 • USA.

Phone (210) 561-9515 • Fax (210) 561-9544

Toll Free (800) 786-5777

Email: service@4adi.com
www.4adi.com

ELISA Kit Components	Amount	Part
Papain Coated Microwell Strip Plate	8-well strips (12)	800-161
Papain Standard, lyophilized	3 vials	800-162
Anti-Papain HRP Conjugate (100X)	0.15 ml	800-164
Sample Diluent (20X)	10 ml	SD20T
Wash Solution Concentrate (100X)	10 ml	WB-100
TMB Substrate	12 ml	80091
Stop Solution	12 ml	80101
Product Manual	1 ea	M-800-160-CPP