

ELISA kits available from ADI (see details at the web site)

- | | | | |
|---------------------|---|--------------|---------------------------|
| #0010 | Human Leptin | | |
| #200-120-AGH | Human globular Adiponectin (gAcrp30) | | |
| #0700 | Human Sex Hormone Binding Glob (SHBG) | | |
| #0900 | Human IGF-Binding Protein 1 (IGFBP1) | | |
| #1000 | Human C-Reactive Protein (CRP) | | |
| #100-110-RSH | Human Resistin /FIZZ3 | | |
| #100-140-ADH | Human Adiponectin (Acrp30) | | |
| #100-160-ANH | Human Angiogenin | | |
| #100-180-APH | Human Angiopoietin-2 (Ang-2) | | |
| #100-190-B7H | Human Bone Morphogenic Protein 7 (BMP-7) | | |
| #1190 | Human Serum Albumin | #1200 | Human Albumin (Urinary) |
| #1750 | Human IgG (total) | #1760 | Human IgM |
| #1800 | Human IgE | #1810 | Human Ferritin |
| #1210 | Human Transferrin (Tf) | #0020 | Beta-2 microglobulin |
| #1600 | Human Growth Hormone (GH) | | |
| #0060 | Human Pancreatic Colorectal cancer (CA-242) | | |
| #1820 | Human Ovarian Cancer (CA125) | #1830 | Human CA153 |
| #1840 | Human Pancreatic & GI Cancer (CA199) | | |
| #1310 | Human Pancreatic Lipase | | |
| #1400 | Human Prostatic Acid Phosphatase (PAP) | | |
| #1500 | Human Prostate Specific Antigen (PSA) | #1510 | free PSA (fPSA) |
| #0500 | Human Alpha Fetoprotein (AFP) | | |
| #0050 | Human Neuron Specific Enolase (NSE) | | |
| #0030 | Human Insulin | #0040 | Human C-peptide |
| #0100 | Human Luteinizing Hormone (LH) | | |
| #0200 | Human Follicle Stimulating Hormone (FSH) | | |
| #0300 | Human Prolactin (PRL) | | |
| #0400 | Human Chorionic Gonadotropin (HCG) | #0410 | HCG-free beta |
| #0600 | Human Thyroid Stimulating Hormone (TSH) | | |
| #1100 | Human Total Thyroxine (T4) | #1110 | Human Free T4 (fT4) |
| #1650 | Human free triiodothyronine (fT3) | #1700 | Human T3 (total) |
| #1850 | Human Cortisol | #1860 | Human Progesterone |
| #1865 | Human Pregnenolone | #1875 | Human Aldosterone |
| #1880 | Human Testosterone | #1885 | Human free Testosterone |
| #1910 | Human Androstenedione | #1920 | Human Estradiol |
| #1925 | Human Estrone | #1940 | Dihydrotestosterone (DHT) |
| #1950 | Human DHEA-sulphate (DHEA-S) | | |
| #3400 | Human serum Neopterin | | |
| #3000 | Human Rheumatoid Factors IgM (RF) | | |
| #3100 | Human anti-dsDNA | | |
| #3200 | Anti-Nuclear Antibodies (ANA) | | |

Alpha Diagnostic Intl. Inc. www.4adi.com 1100/151113AA

TOTAL THYROXINE (TOTAL T₄)

ELISA KIT Cat. No. 1100, 96 Tests

For Quantitative Determination of Total T₄ In Human Serum



For In Vitro Research Use Only



4638 N Loop 1604 West • San Antonio • Texas 78249 • USA.
Phone (210) 561-9515 • Fax (210) 561-9544
Toll Free (800) 786-5777

Email: Techsupport@4adi.com

Web Site: www.4adi.com

ELISA KIT Cat. No. 1100 (96 tests)

For Quantitative Determination of Total T₄ in Human Serum

Kit Components (96 tests)	Cat #
streptavidin coated Strip plate, (96 wells)	1101
T4 Standard A , 0.5 ml; 0 ug/dL	1102A
T4 Standard B , 0.5 ml; 2.0 ug/dL	1102B
T4 Standard C , 0.5 ml; 5.0 ug/dL	1102C
T4 Standard D , 0.5 ml; 10.0 ug/dL	1102D
T4 Standard E , 0.5 ml; 15.0 ug/dL	1102E
T4 Standard F , 0.5 ml; 25.0 ug/dL	1102F
Assay Diluent , 7 ml	1105AD
T4- HRP Conjugate conc. (11X), 0.7 ml,	1103
Anti-T4 Biotin Solution ; 7 ml	1104
HRP Substrate Solution, 12 ml	1100-SS
Wash Buffer (20X) ; 25 ml, dilute with 475 ml of distilled water	1 1 0 0 - W B
Stop Solution , 12 ml	1100-ST
Complete Instruction Manual	M-1100

Intended Use:

The Thyroxine (T₄) ELISA Kit is intended for the detection of Total T₄ in human serum or plasma. For **In Vitro research use only (RUO)**.

Introduction

The thyroid gland produces T₄, triiodothyronine T₃ and calcitonin. The first two hormones are synthesized by the gland following entrapment of iodine, conversion to iodine, and coupling of iodine with tyrosine, followed by coupling of two iodinated tyrosine molecules. T₄ and T₃ so formed are attached to thyroglobulin for storage and are released, as needed, as protease splits them from the globulin. Thyroxine is a highly active thyrometabolic hormone, which exists, in protein-bound and unbound forms. For T₄ can be measured more easily and with greater accuracy than T₃, determination of total T₄ by immunoassay is the most reliable and convenient screening test available for detecting thyroid disorders in man.

Release of T₄ and T₃ from the thyroid is greatly influenced by pituitary-thyroid stimulating hormone (TSH), which in turn is influenced by hypothalamic thyrotropin-releasing hormone (TRH). Normally, increased blood levels of T₄ and T₃ act to decrease the amount of TSH secreted, thereby reducing the production and release of T₄ and T₃. Decreased blood levels of T₄ and T₃ produce the opposite effect, leading to increased production and secretion of T₄ and T₃. In this manner a normal circulating thyroid hormone balance is maintained, circulating T₄ and T₃ are bound largely to thyroxine binding globulin (TBG). To a lesser extent they are bound to thyroxine binding prealbumin (TBPA) and, when present in excess, to albumin. Usually T₄ and T₃ concentration ratio is about 9:1, however T₃ has considerably greater physiological activity. It is the small free fraction (0.1% of the total or less) that is physiologically active and determines the clinical thyroid status of the patient's hyperthyroid, euthyroid, or hypothyroid.

Inter-assay precision:

Three serum samples (4-15 ug/dL) were run in duplicate in sixteen independent assays. The samples showed good inter-assay precision (2.7-4 % CV).

3. RECOVERY

A known amount of total T₄ (3-24 ug/dL) was added to three patient sera (with original total T₄ concentrations of 3 and 6 ug/dL) and the final total T₄ concentration measured. The assay showed good mean recoveries of about 93% (range 85-96%).

4. LINEARITY

A serum sample containing 24 ug/dL was diluted with a series of T₄-free serum. The dilutions were tested and the T₄ recoveries were compared with the expected concentration. The samples showed excellent mean recoveries of about 102% (range 90-109%).

5. SPECIFICITY

The specificity of total T₄ ELISA kit was determined by measuring interference from high concentrations of 3,5-diiodothyronine (up to 10 ug/dL), 3,3',5-triiodothyronine (rT₃, up to 1000 ng/ml), 3,3',5-triiodothyronine (T₃, up to 59 nmol/l), 3,3',5-triiodothyroacetic acid (up to 42 nmol/l), 3,3',5-triiodothyropropionic acid (up to 63 nmol/l), Aspirin (10 mg/dL) iodoacetic acid (10 ug/dL), phenylbutazone (10 mg/dL). No significant cross-reaction was observed with any of these compounds.

6. CORRELATIVE STUDY

ADI total T₄ elisa kit was compared with Inc Star Clinical assay T₄ RIA by analyzing 99 samples (2.6-25.5 ug/dL). The regression analyses showed good correlation between the two assays.

	Inc Star RIA	ADI T4 ELISA
Mean Sample (ug/dL)	8.682	8.854
Slope = 0.98288; Intercept =0.790		Correlation Coefficient =0.852
ADI = 0.929 (Inc stat T4)+ 0.790		

7. SPECIES REACTIVITY AND PUBLICATIONS

T₄ is the same in human and animals. ADI human T₄ elisa Kit has been used in mouse and rat samples (see below).


Maglich 2004, J. Biol. Chem., 279, 19832 – 19838, mice serum T₄
Huang W, 2005, J. Nutr., 135: 1631 - 1635 T₄ ELISA
Xiao CW, 2004, J. Nutr., 134: 743 – 749, T₄ in rat Plasma

General References: Ekins R et al (1969) Clin. Biochem. 2; Skelly D et al (1973) Clin. Chem. 2, 19; Shcuurs A et al (1977) Clin. Chim. Acta 5, 81; Scharpe S et al (1976) Clin. Chem. 6, 22; Wisdom G et al (1976) Clin Chem. 8, 22, Schall R et al (1978) Clin. Chem. 10, 24;

WORKSHEET OF TYPICAL ASSAY

Wells	Stds/samples	Mean A ₄₅₀ nm	A/A0 x 100 Index	Calculated Concn (ug d/L)
A1, A2	Std. A (0 ug/dL)	2.679		
B1, B2	Std. B (2 ug/dL)	1.961		
C1, C2	Std. C (5 ug/dL)	1.322		
D1, D2	Std. D (10 ug/dL)	0.901		
E1, E2	Std. E (15 ug/dL)	0.659		
F1, F2	Std. F (25 ug/dL)	0.396		
G1, G2	Sample 1			

NOTE: These data are for demonstration purpose only. A complete standard curve must be run in every assay to determine sample values. Each laboratory should determine their own normal reference values.



A typical std. assay curve (do not use this for calculating sample values)

PERFORMANCE CHARACTERISTICS

1. DETECTION LIMIT

Based on sixteen replicates determinations of the zero standard, the minimum concentration of total T₄ detected using this assay is 0.5 ug/dL. The detection limit is defined as the value deviating by 2 SD from the zero standard.

2. PRECISION

Intra-assay precision:

Three serum samples (mean total T₄ concentrations 3.5, 6.8, 13.9 ug/dL) were run in five separate runs. The samples showed good intra-assay precision with %CV of 9-11.

Alpha Diagnostic Intl. Inc. www.4adi.com 1100/151113AA Page 5

PRINCIPLE OF THE TEST

Total T₄ ELISA kit is based on competitive binding of human thyroxine from serum samples and enzyme-labeled T₄ to T₄-specific antibodies immobilized on microtiter well plates. In the assay, total T₄ is released from its binding proteins by a releasing agent present in the assay buffer. After a washing step, chromogenic substrate is added and color developed. The enzymatic reaction (blue color) is inversely proportional to the amount of T₄ present in the sample. The reaction is terminated by adding stopping solution (converts blue to yellow). Absorbance is then measured on a microtiter well ELISA reader at 450 nm. and the concentration of T₄ in samples and control is read off the standard curve.

MATERIALS AND EQUIPMENT REQUIRED

Adjustable micropipet (20-100 ul) and multichannel pipet with disposable plastic tips. Reagent troughs, plate washer (recommended) and ELISA plate Reader.

PRECAUTIONS

The Alpha Diagnostic International Total T₄ ELISA test is intended for *in vitro research* use only. The reagents contain proclin-300 as preservative; necessary care should be taken when disposing solutions. The Control Serum has been prepared from human sera shown to be negative for HBsAg and HIV antibodies. Nevertheless, such tests are unable to prove the complete absence of viruses; therefore, sera should be handled with appropriate precautions.

Applicable **MSDS**, if not already on file, for the following reagents can be obtained from ADI or the web site.

TMB (substrate), H₂SO₄ (stop solution), and Proclin-300 (0.1% v/v in standards, sample diluent and HRP-conjugates).

All waste material should be properly disinfected before disposal. Avoid contact with the stop solution (1N sulfuric acid).

SPECIMEN COLLECTION AND HANDLING

Collect blood by venipuncture, allow to clot, and separate the serum by centrifugation at room temperature. Do not heat inactivate the serum.. If sera cannot be immediately assayed, these could be stored at -20°C for up to six months. Avoid repeated freezing and thawing of samples. No preservatives should be added to the serum.

REAGENTS PREPARATION

Dilute wash buffer (1:20) with distilled water (25 ml stock in 475 ml). Store at 4oC

Enzyme conjugate (11X) dilute 1:11 with assay diluent. (For example, dilute 160µl of enzyme conjugate with 1.6ml of assay diluent for 16 wells (A slight excess of solution is made). This reagent should be used within twenty-four hours for maximum performance of the assay. Store at 2-8°C. **Dilute conjugate in required amounts only and do not store diluted conjugate.**

Alpha Diagnostic Intl. Inc. www.4adi.com 1100/151113AA Page 2

STORAGE AND STABILITY

The microtiter well plate and all other reagents are stable at 2-8°C until the expiration date printed on the label. The whole kit stability is usually 6 months from the date of shipping under appropriate storage conditions. The unused portions of the standards should be stored at 2-8°C or stored frozen in small aliquots.

TEST PROCEDURE (ALLOW ALL REAGENTS TO REACH ROOM TEMPERATURE BEFORE USE).

Remove required number of coated strips and arrange them on the plate. Store unused strips in the bag. Dilute conjugate in required amounts. **Dilute wash buffer (1:20) with distilled water (25 ml stock in 475 ml).** Store at 4°C.

1. Label or mark the microtiter well strips to be used on the plate. Store unused strips in the sealed pouch at 4°C.
2. Pipet **25 ul of standards**, control, and serum samples into appropriate wells in *duplicate*.
3. Add **50 ul of enzyme conjugate** into each well. Add **50 ul of T4-Antibody biotin Solution** to all wells. Mix gently for 20-30 seconds. Cover the plate and incubate for **60 minutes** at room temperature.
4. Aspirate and **wash the wells 3 times** with **300 ul** of wash buffer. We recommend using an automated ELISA plate Washer for better consistency. Failure to wash the wells properly will lead to high blank or zero values. If washing manually, plate must be tapped over paper towel between washings to ensure proper washing.
5. Dispense **100 ul TMB substrate per well**. Mix gently for 5-10 seconds.
6. Cover the plate and incubate for **15 minutes** at room temperature. Blue color develops in standards and samples.
7. Stop the reaction by adding **50 ul of stop solution** to all wells at the same timed intervals as in step 6. Mix gently for 5-10 seconds. Blue color turns yellow.
8. Measure the **absorbance at 450 nm** using an ELISA reader within 30 min.

NOTES

Read instructions carefully before the assay. Do not allow reagents to dry on the wells. Careful aspiration of the washing solution is essential for good assay precision. Since timing of the incubation steps is important to the performance of the assay, pipet the samples without interruption and it should not exceed 5 minutes to avoid assay drift. If more than one plate is being used in one run, it is recommended to include a standard curve on each plate. The unused strips should be stored in a sealed bag at 4°C.

Addition of the HRP substrate solution starts a kinetic reaction, which is terminated by dispensing the stopping solution. Therefore, keep the incubation time for each well the same by adding the reagents in identical sequence. Plate readers measure absorbance vertically. Do not touch the bottom of the wells.

Limitations

Serum samples containing more than 24 ug/dL T₄ should be diluted with the zero standard (standard A or normal saline) and the results obtained should be multiplied by the appropriate dilution factor. Samples containing T₄ <0.5 ug/dL are analyzed by diluting with 2 ug/dL T₄ calibrator to extend the curve.

Calibrators and controls from other manufacture should not be used as they may contain serum preservatives incompatible with ADI's ELISA reagents.

Whenever laboratory data conflicts with clinical findings or impressions, clinical judgment should be exercised and additional evaluations undertaken.

Use of ADI's reagent in a study of euthyroid patients in one geographic location will yield a normal range. It is recommended that laboratories adjust normal values to reflect geographic and population differences specific to the patients they serve.

CALCULATION OF RESULTS

Draw the standard curve on a linear graph paper by plotting net absorbance values of standards against appropriate total T₄ concentrations. Read off the total T₄ concentrations of the control and patient samples.

EXPECTED VALUES

In a study of 75 euthyroid patients yielded a normal range of (4-12 ug/dL) at the 95% confidence limit. It is recommended that researchers adjust normal values to reflect geographic and population differences.