

## INTENDED USE

The **Human Anti-HPV11L1 IgG ELISA Kit** is an immunoassay suitable for detecting and quantifying IgG antibody activity specific for **HPV subtype 11 L1** protein in serum or plasma. Other biological fluids, including tissue culture medium, may be validated for use.

This immunoassay is suitable for:

- Determining **immune status** relative to non-immune controls;
- Assessing efficacy of **vaccines**, including dosage, adjuvantcy, route of immunization and timing;
- Qualifying and standardizing vaccine batches & protocols

The assay is for research use only (RUO) and is not intended nor validated for diagnosing HPV. Reagents contain no virus.


## GENERAL INFORMATION

**Human papillomavirus (HPV)** is a group of more than 200 related viruses, of which more than 40 are spread through direct sexual contact. Several HPV types cause genital warts, and about a dozen HPV types can cause certain cancer types – cervical, anal, oropharyngeal, penile, vulvar and vaginal.

The HPV virus is small and non-enveloped, with a single, circular dsDNA contained within a T=7 icosahedral capsid composed of 72 pentamers. The HPV genome is composed of six early (**E1, E2, E3, E4, E6 and E7**) and two late (**L1 and L2**) proteins. The late proteins, translated in the upper layers of the host epithelium, serve as structural proteins that encapsulate the amplified viral genomes.

FDA-approved vaccines that have shown to produce protective antibodies include **Gardasil** (types 6,11,16,18) & **Gardasil 9** (with add'l types 31,33,45,52 & 58) and **Cervarix** (types 16,18), which are composed of recombinant L1 protein as virus-like particles (VLP).

## PRINCIPLE OF THE TEST



The Human Anti-HPV11L1 IgG ELISA kit is based on the binding of antibody in samples to HPV11L1 (serotype specific) antigen immobilized on the microwells, and bound antibody is detected by antibody-HRP conjugate. After a washing step, chromogenic substrate (TMB) is added and color is developed by the enzymatic reaction of HRP on

the substrate, which is directly proportional to the amount of antibody present in the sample. Stopping Solution is added to terminate the reaction, and absorbance at 450nm is then measured using an ELISA microwell reader. The presence of Human HPV11L1-IgG antibody in samples is determined relative to anti-HPV11L1 Calibrators supplied in the kit.

## PRODUCT SPECIFICATIONS

### Specificity

Purified recombinant HPV11L1 protein virus-like particles (VLP) (baculovirus-insect cells; 496-aa, ~55 kDa) are used to coat the microwells; thus, no other antibody specificity is detectable in the assay. HPV11L1 shares substantial sequence homology with other HPV11 subtypes. The Anti-Human IgG HRP conjugate reacts specifically with human IgG class antibodies; IgA, IgM and IgE antibody would not be measured above background signals.

## KIT CONTENTS

The microtiter well plate and all other reagents, if unopened, are stable at 2-8°C until the expiration date printed on the box label. Stabilities of the working solutions are indicated under Reagent Preparation.

**To Be Reconstituted:** Store as indicated.

Component	Preparation Instructions
<b>Wash Solution Concentrate (100x)</b> Cat. No. WB-100, 10ml	Dilute the entire volume 10ml + 990ml with distilled or deionized water into a clean stock bottle. Label as <b>Working Wash Solution</b> and store at 4°C for long term and ambient temp. for short term.
<b>Sample Diluent Concentrate (20x)</b> Cat. No. SD-20T, 10ml	Dilute the entire volume, 10ml + 190ml with distilled or deionized water into a clean stock bottle. Label as <b>Working Sample Diluent</b> and store at 2-8°C until the kit lot expires or is used up.
<b>Anti-Human IgG-HRP Conjugate Concentrate (100x)</b> Part: H-HuG.2a11, 0.15ml	Peroxidase conjugated anti-human IgG in buffer with detergents and antimicrobial as stabilizers. Dilute fresh as needed; 10ul of concentrate to 1ml of <b>Working Sample Diluent</b> is sufficient for 1 8-well strip. Use within the working day and discard. Return 100X to 2-8°C storage.

**Ready For Use:** Store as indicated on labels.

Component	Part	Amt	Contents
<b>HPV11 Antigen Coated Strip Plate</b>	550-HPV11	8-well strips (12)	Coated with recombinant HPV11L1 protein, and post-coated with stabilizers.
<b>Anti-HPV11 Calibrators</b>			
1 U/ml	550212B	0.65 ml	Four (4) vials, each containing anti-HPV11L1; in buffer with antimicrobial as stabilizers.
2.5 U/ml	550212C	0.65 ml	
5 U/ml	550212D	0.65 ml	
10 U/ml	550212E	0.65 ml	
<b>Anti-HPV11 Positive Control</b>	550-212PC	0.65 ml	Serum with anti-HPV11L1 reactivity; [Value Range on Label]
<b>Low NSB Sample Diluent</b>	TBTm <b>Not for HRP dilution</b>	30 ml	Buffer with protein, detergents and antimicrobial. Use as is for sample dilution
<b>TMB Substrate</b>	80091	12 ml	Chromogenic substrate for HRP containing TMB and peroxide.
<b>Stop Solution</b>	80101	12 ml	Dilute sulfuric acid.

### Materials Required But Not Provided:

- Pipettors and pipettes that deliver 100ul and 1-10ml.
- Disposable glass or plastic 5-15ml tubes for diluting samples and Anti-Human IgG HRP Concentrate.
- Stock bottle to store diluted Wash Solution; 0.2 to 1L.
- Distilled or deionized water to dilute reagent concentrates.
- Microwell plate reader at 450 nm wavelength and ELISA plate washer

## ASSAY DESIGN AND SET-UP

### Sample Collection and Handling

Serum and other biological fluids may be used as samples with proper dilution to avoid solution matrix interference. For **serum**, collect blood by venipuncture, allow clotting, and separate the serum by centrifugation at room temperature. For other samples, clarify the sample by centrifugation and/or filtration prior to dilution in Sample Diluent. If samples will not be assayed immediately, store refrigerated for up to a few weeks, or frozen for long-term storage.

### Antibody Stability & Dilution

Initial dilution of serum into **Working Sample Diluent** (WSD) is recommended to stabilize antibody activity. This enhances reproducible sampling, and stabilizes the antibody activity for years, stored refrigerated or frozen. Further dilution into **Low NSB Sample Diluent** (LNSD), which provides the lowest assay background, should be at least 10 times the initial dilution and performed the same day as the assay.

Example: Initial (1:5): **10ul** serum + **40ul** WSD

Further (1:50): **10ul** initial (1:5) + **90ul** LNSD (1:50)

**Note:** Samples may be diluted straight into the Low NSB Sample Diluent without initial dilution into WSD; stable refrigerated for a month.

### Assay Design

Review Interpretation of Results and Limits of the Assay (p5-7) before proceeding:

- Select the proper sample dilutions accounting for expected potency of positives and minimizing non-specific binding (NSB) and other matrix effects; for example, net signal for non-immune samples should be lower than the **1 U/ml Calibrator**. This is usually 1:200 or greater dilution for human serum/plasma with normal levels of IgG and IgM.
- Run the **Anti-HPV11 Positive Control**; value range on label.
- Run a Sample Diluent **Blank**. This signal is an indicator of proper assay performance, especially of washing efficacy, and is used for net OD calculations, if required. Blank OD should be <0.3.
- Run a set of **Calibrators**, which validate that the assay was performed to specifications: **10 U/ml** should give a high signal (>1.5 OD); **1 U/ml** should give a low signal which can be used to discriminate at the Positive/Negative threshold (see Interpretation of Results, p. 5).

### Plate Set-up

Bring all reagents to room temperature (18-30° C) equilibration (at least 30 minutes).

- Determine the number of wells for the assay run. Duplicates are recommended, including 8 Calibrator wells and 2 wells for each sample control to be assayed.
- Remove the appropriate number of microwell strips from the pouch and return unused strips to the pouch. Reseal the pouch and store refrigerated.
- Add 200-300ul Working Wash Solution to each well and let stand for about 5 minutes. Aspirate or dump the liquid and pat dry on a paper towel before sample addition.

## Assay Procedure

ALL STEPS ARE PERFORMED AT ROOM TEMPERATURE. After each reagent addition, gently tap the plate to mix the well contents prior to beginning incubation.

### 1. 1<sup>st</sup> Incubation [100ul – 60 min; 4 washes]

- Add 100ul of calibrators, samples and controls each to pre-determined wells.
- Tap the plate gently to mix reagents and incubate for 60 minutes.
- Wash wells 4 times and pat dry on fresh paper towels. As an alternative, an automatic plate washer may be used. Improper washes may lead to falsely elevated signals and poor reproducibility.

### 2. 2<sup>nd</sup> Incubation [100ul – 30 min; 5 washes]

- Add 100ul of diluted Anti-Human IgG HRP to each well.
- Incubate for 30 minutes.
- Wash wells 5 times as in step 2.

### 3. Substrate Incubation [100ul – 15 min]

- Add 100ul TMB Substrate to each well. The liquid in the wells will begin to turn blue.
- Incubate for 15 minutes in the dark, e.g., place in a drawer or closet.

Note: If your microplate reader does not register optical density (OD) above 2.0, incubate for less time, or read OD at 405-410 nm (results are valid).

### 4. Stop Step [Stop: 100ul]

- Add 100ul of Stop Solution to each well.
- Tap gently to mix. The enzyme reaction will stop; liquid in the wells will turn yellow.

### 5. Absorbance Reading


- Use any commercially available microplate reader capable of reading at 450nm wavelength. Use a program suitable for obtaining OD readings, and data calculations if available.
- Read absorbance of the entire plate at 450nm using a single wavelength within 30 minutes after Stop Solution addition. If available, program to subtract OD at 630nm to normalize well background.

## INTERPRETATION OF RESULTS

### Method A. Antibody Activity Threshold Index

Compare Samples to 1 U/ml Calibrator or Internal Control  
= Positive/Negative Cut-off.

#### Example:



#### Results

The **sensitivity** of the assay to detect anti-HPV11L1 IgG, from either natural exposure or vaccination, is controlled so that the 1 U/ml Calibrator represents a threshold OD for most true positives in human serum diluted to 1:200 or less. Visual inspection of the data in the above graph shows the following:

**Calibrators** – dilution curve of an anti-HPV antibody, derived from HPV11L1 vaccination, shows the OD range of the assay; high value indicates optimal sensitivity of the assay.

**1 U/ml:** a 'Cut-off' line has been drawn to indicate a threshold distinguishing between **Positive/Negative**. This is not a clear-cut threshold, rather a low OD area that could represent either low positives or high background negatives.

**Positive Control** – an anti-HPV11L1 serum; value range is on the label. The control can be used to gauge reproducibility and to normalize between-assay variation.

**Internal Control** – a low level positive from an immunized individual that represents the lab's experience in distinguishing low positive from negative samples. This should be run in each assay to supplement the 1 U/ml Calibrator for Positive/Negative discrimination purposes.

**Samples A,B,C,D** – 3 samples (A, B, C) are negative: below the threshold; 1 sample (D) is positive: clearly above the threshold.

The 1 U/ml Calibrator can be used to calculate a **Threshold Index** that numerically discriminates Positive/Negative:


- ❖ Divide each Sample net OD by the 1 U/ml Calibrator net OD. Values above 1.0 are a measure of **Positive** Antibody Activity; below 1.0 are **Negative** for antibody.

## INTERPRETATION OF RESULTS (cont)

### Example:

#### Human Serum IgG

A panel of human sera of unknown history was tested for anti-HPV11 IgG (1:200 dilution in Low NSB Sample Diluent). **Threshold Index** was calculated using the 1 U/ml Calibrator.



#### Results

**Anti-HPV IgG:** 9 of 10 samples were negative at 1:200 dilution (below the 1.0 Threshold Index); one sample (J) was borderline.

#### Notes:

- Positives** may be due to prior encounter with HPV, from exposure to an antigen with common epitopes, or from immunization.
- When the **Positive Index** is **above 5.0**, using a dilution curve to calculate titer is a more accurate quantitation method (see Method C).
- The **sensitivity** of the assay may be adjusted by changing the sample dilutions: a) increase dilution (e.g., 1:200) to lower the signals of borderline positives to negative; b) decrease dilution (e.g., 1:50) to convert borderline samples to positive. With the latter, the values of negatives may increase, so an alternative threshold should be considered using known negatives to develop a **Positive Index** (see below) or use an **Internal Control** (Page 5).

#### B. Positive Index


Experimental sample values may be expressed relative to the values of Control or Non-immune samples, by calculation of a **Positive Index**. One typical method is as follows:

- Calculate the net OD mean + 2 SD of the Control/Non-immune samples = **Positive Index**.
- Divide each sample net OD by the Positive Index. Values above 1.0 are a measure of **Positive** Antibody Activity; below 1.0 are **Negative** for antibody.

## INTERPRETATION OF RESULTS (cont)

### C. Antibody Titer

The most accurate method for comparing antibody potencies is by calculation of a titer, using an OD reading midrange in the dilution curves of each antibody as **Index**. In the example below, **IgG** titers were calculated as inverse of the dilution that produced a **1.0 OD** in the assay.



#### Results

Antisera are from individuals immunized with Gardasil (Merck) HPV 6+11+16+18.

**Mouse Anti-HPV:** Pool of 5 hyperimmunized mice; Titer = **130k**.

**Rabbit Anti-HPV:** Hyperimmunized individual; Titer = **1.4 kk**.

**Human Anti-HPV:** WHO International Standard for HPV16; Titer = **375 = 13.3 mIU/ml**.

#### Limits of the Assay

- The assay detects and quantifies IgG antibodies directed to the L1 protein. Patients may have HPV infection without producing antibodies specific to L1.
- Anti-HPV antibody levels of an infected patient may be below detection threshold related to the time course of the infection, e.g., too early for positive titer development.
- Samples from non-patients may be elevated due to prior exposure to the human papilloma virus.

## PRECAUTIONS AND SAFETY INSTRUCTIONS

Calibrators, Sample Diluent, and Antibody HRP contain bromonitrodioxane (BND: 0.05%, w/v). Stop Solution contains dilute sulfuric acid. Follow good laboratory practices, and avoid ingestion or contact of any reagent with skin, eyes or mucous membranes. All reagents may be disposed of down a drain with copious amounts of water. MSDS for TMB, sulfuric acid and BND can be requested or obtained from the ADI website: [http://4adi.com/commerce/info/showpage.jsp?page\\_id=1060&category\\_id=2430&visit=10](http://4adi.com/commerce/info/showpage.jsp?page_id=1060&category_id=2430&visit=10)

# Human Anti-HPV11L1 IgG ELISA Kit

**Cat. # 550-111-PHG, 96 tests**

**For the Detection and Quantitation of Anti-HPV11L1 (Late Protein: type 11) IgG in Human Serum/Plasma**

*For research use only, not for diagnostic or therapeutic use.*



**ALPHA DIAGNOSTIC INTERNATIONAL**

4638 N Loop 1604 West Center Drive • San Antonio • Texas 78249 •  
USA, Phone (210) 561-9515 • Fax (210) 561-9544  
Toll Free (800) 786-5777 Email: [service@4adi.com](mailto:service@4adi.com)

ELISA Kit Components	Amount	Part
HPV11 Antigen	8-well strips (12)	550-HPV11
Anti-HPV11 Positive Control	0.65 ml	550-212PC
Anti-HPV11 Calibrator 1 U/ml	0.65 ml	550-212B
Anti-HPV11 Calibrator 2.5 U/ml	0.65 ml	550-212C
Anti-HPV11 Calibrator 5 U/ml	0.65 ml	550-212D
Anti-HPV11 Calibrator 10 U/ml	0.65 ml	550-212E
Anti-Human IgG HRP Conjugate (100X)	0.15 ml	H-HuG.2a11
Sample Diluent (20x)	10 ml	SD20T
Low NSB Sample Diluent	30 ml	TBTm
Wash Solution Concentrate (100X)	10 ml	WB-100
TMB Substrate	12 ml	80091
Stop Solution	12 ml	80101
Product Manual	1 ea	550-111-PHG